


Screen all your COPD patients with

AlphaID™



Cheek Swab



Dried Blood Spot

AlphaID provides *completely free* screening for
alpha₁-antitrypsin deficiency (alpha-1)

Alpha-1 is the #1 known genetic risk factor for COPD¹

GRIFOLS

Alpha-1 is the #1 Known Genetic Risk Factor for COPD¹



Alpha-1 may be a contributing cause in **up to 3% of COPD cases** in the United States—**up to 300,000 people**¹



First described in 1963, alpha-1 is **caused by a mutation** of the SERPINA1 gene located on chromosome 14²




Alpha₁ antitrypsin protects the lungs from proteolytic attack, particularly by neutrophil elastase. Without sufficient levels of alpha₁ antitrypsin, lung damage occurs¹



A landmark study indicated that **severe alpha₁ significantly shortened life span**²

- Mean survival of the alpha-1 cohort (n=120) was 16% at 60 years of age, compared with 85% for the general population



Study Design, 1988: All 120 Pi ZZ subjects evaluated as inpatients in NIH Clinical Center. Various demographic features, clinical findings, and electrocardiogram evaluations were obtained from medical records. Blood and urine laboratory findings were obtained through computerized archives of the Clinical Center Laboratories. Chest radiographs and scintigraphic scans were evaluated de novo at the time of this analysis, with knowledge of the diagnosis but without knowledge of any of the other data. Lung function tests were obtained from the computerized files of the Pulmonary Branch. Mortality data were obtained at the time of this analysis. Alpha₁ antitrypsin protein phenotypes were determined by a combination of isoelectric focusing of serum, quantitation of AAT levels in serum, and family studies.²

*Reprinted with permission of the American Thoracic Society. Copyright 2012, American Thoracic Society.

Alpha-1 is Underdiagnosed¹



More than 90% of people with alpha-1 are estimated to be undiagnosed




The **average delay** between the onset of symptoms and diagnosis of alpha-1 is **8.3 years**



The average **number of physicians** seen before being correctly diagnosed is **2.7**

Alpha-1 occurs in both nonsmokers and smokers



*378 Pi ZZ patients.

[†]Stage III severe COPD, GOLD guidelines.

These Major Respiratory Organizations Recommend Targeted Genetic Screening for Alpha-1 in All Patients With COPD⁵⁻⁹



ATS: American Thoracic Society



CHEST: American College of Chest Physicians



COPD Foundation: Chronic Obstructive Pulmonary Disease

GOLD

GOLD: Global Initiative for Chronic Obstructive Lung Disease

WHO

WHO: World Health Organization

The American Thoracic Society recommends screening for alpha-1 regardless of age, smoking history, or FEV₁ in⁵:

- All adults with symptomatic COPD or asthma
- All adults with symptomatic asthma whose airflow obstruction is incompletely reversible after bronchodilator therapy
- Asymptomatic patients with persistent obstruction on pulmonary function tests and with identifiable risk factors (eg, smoking, occupational exposure)
- All immediate family members of patients who are diagnosed with alpha-1



Screening patients with COPD for alpha-1:

- Identifies those with a severe alpha-1 for whom treatment may be appropriate⁵
- Identifies carriers of deficient alleles, which can be passed onto children¹
- Provides incentive for smoking cessation¹⁰
 - In a follow-up study, 59% of patients with severe alpha-1 attempted to quit smoking after receiving test results and minimal counseling

Use AlphaID to Screen for Alpha-1



Each sample is processed by Matrix Clinical Labs, a **CLIA-certified**, independent laboratory and **analyzed by alpha-1 experts**¹¹



Comprehensive results are available within **5 business days** from the time the sample is returned to the laboratory¹¹

- Results are accessible through the confidential, easy-to-use portal at myalphaID.com
- Results are also received by fax



Alpha-1 expert clinicians are available for **personal consultation to discuss abnormal results**



Results are **strictly confidential** and are shared only with the ordering clinician and designated office staff



Provided by Grifols at **no cost** and may not be billed to patients or their insurers



Cheek Swab



Dried Blood Spot



A free, fast, and easy way to rule out alpha-1



	AlphaID:
Fast	A gentle cheek swab that only takes seconds. There are no chemicals to add and no cards to dry
Familiar	Uses the same DNA collection and screening methods as many at-home genetic testing kits
Reliable	Delivers accurate and reliable genotyping results ¹²
Specific	Screens for the 14 most-prevalent mutations associated with alpha-1, including S, Z, F, I, rare, and null alleles ¹²
Confidential	Results are available only to the ordering clinician and designated office staff
Completely Free	Provided by Grifols and may not be billed to patients or their insurers

The AlphaID kit contains everything you need to collect and submit a cheek swab sample for alpha-1 screening

- Instructions
- Requisition form
- Cheek swab
- Specimen bag
- Mailer box
- Prepaid shipping label



EXAMPLE RESULTS REPORT

Patient Report

Test, Grifols Z/Z
 DOB: 1/1/2000
 MRN: 23-033-0000778
 GENDER: M

Ordering Provider: **Test, Grifols**
 Collected: 3/3/2023
 Received: 3/3/2023
 Reported: 3/3/2023

AlphaID

Specimen ID: T2306201131

ABNORMAL TEST RESULT

Test	Result	Units	Flag	Reference
Alpha-1 Antitrypsin Genotype Test	Z/Z	N/A	A	M/M

Interpretation:

PI^{*}M/M is the normal genotype. Abnormal genotype results are associated with increased risk of alpha-1 antitrypsin deficiency. Clinical correlation with patient's serum alpha-1 antitrypsin level and clinical presentations are recommended.

The alpha-1 antitrypsin genotype test specifically detects 14 common mutations associated with alpha-1 antitrypsin deficiency. No other mutations are detected. This test is FDA cleared and its performance characteristics determined by Trillium Health.

The laboratory can be reached at alphaid@trillium-health.com or (888)261-2812, option 2.

AlphaID reports the alpha-1 genotype. Any result other than a homozygous MM is considered abnormal.



Confirming an alpha-1 diagnosis is simple



	AlphaID CONFIRM:
Efficient	Uses a dried blood spot card to screen for alpha-1 or confirm the results from an AlphaID cheek swab
Comprehensive	Confirms alpha-1 genotype and reports the serum level of alpha ₁ antitrypsin
Reliable	Delivers accurate and reliable genotyping results ¹²
Specific	Screens for the 14 most-prevalent mutations associated with alpha-1, including S, Z, F, I, rare, and null alleles ¹²
Confidential	Results are available only to the ordering clinician and designated office staff
Completely Free	Provided by Grifols and may not be billed to patients or their insurers

The AlphaID CONFIRM kit contains everything you need to collect and submit a dried blood sample for alpha-1 screening

- Instructions
- Requisition form
- Two lancets
- Alcohol pads, gauze, and bandage
- Dried blood spot card
- Postage-paid envelope



EXAMPLE RESULTS REPORT

Patient Report

Test, Grifols_Confirm
 DOB: 2/2/2000
 MRN: 23-079-0000211
 GENDER: F

Ordering Provider: **Test, Grifols**
 Collected: 3/24/2023
 Received: 3/24/2023
 Reported: 3/27/2023

AlphaID Confirm

Specimen ID: T2308303563

ABNORMAL TEST RESULT

Test	Result	Units	Flag	Reference
Alpha-1 Antitrypsin, Quantitation	250	mg/dL		90-200
Alpha-1 Antitrypsin Genotype Test	M/Z	N/A	A	M/M

Interpretation:

M/M is the normal genotype. Abnormal genotype results are associated with increased risk of alpha-1 antitrypsin deficiency. Clinical correlation with patient's serum alpha-1 antitrypsin level and clinical presentations are recommended.


The alpha-1 antitrypsin genotype test specifically detects 14 common mutations associated with alpha-1 antitrypsin deficiency. No other mutations are detected. This test is FDA cleared and its performance characteristics determined by Trillium Health.

AlphaID CONFIRM reports the alpha-1 genotype as well as the alpha₁-antitrypsin level (normal is 90.0 to 200.0 mg/dL).

Abnormal Alpha-1 Alleles Are Common¹¹

Genetic screening found ~1 in 7 patients had deficient alleles

- Of screened patients, 15% were positive for a deficient allele
- F and I allele combinations, as well as rare variants, were found in >1% of >1 million patients screened




AlphaID™ and AlphaID CONFIRM detect 14 genetic variations associated with alpha-1¹¹

- **Most common normal allele:** M*
- **Most common deficient alleles:** S, Z, F, I
- **Less frequent deficient alleles:** P Lowell, P Duarte, QO Cardiff, Y Barcelona, M Procida, M Malton, M Palermo, M Nichinan, S Iiyama, M Heerlen
- **Null alleles:** Q0 Granite Falls, Q0 West, Q0 Bellingham, Q0 Mattawa, Q0 Ourem, Q0 Clayton, Q0 Saarbruecken

*M allele means that none of the 14 allelic variants interrogated by the A1AT Genotyping test are detected in the SERPINA1 gene, but other variants could be present. Notes: The various background M alleles (M1,M2,M3,M4) are not detected by PCR but explain the different nomenclatures used for a same mutation.

All deficient alleles compromise alpha₁ antitrypsin production, but the effect is most pronounced for the Z allele¹¹




The My AlphaID Portal is a Convenient Way to Manage Alpha-1 Results for All of Your Patients

AlphaID kits are processed by a CLIA-certified, independent laboratory

My AlphaID Portal access is granted by following this 3-step process:

1. Visit myalphaid.com and click 'sign up now'.
2. Provide email address and request verification code.
3. Once received, enter verification code, choose your password, and choose account type for your practice.


You will be asked to open an individual account or group account manager. If you have multiple provider NPI numbers to enter, please choose group account manager.



My AlphaID Portal provides:

- Convenient Management
- Comprehensive Reporting
- Timely Results
- Staff Access

Access confidential results processed by a certified, independent laboratory



Use AlphaID to Screen Appropriate Patients for Alpha-1

- AlphaID provides **completely free** screening for alpha₁-antitrypsin deficiency (alpha-1), the #1 known genetic risk factor for COPD¹
- More than 90% of people with alpha-1 are estimated to be undiagnosed¹
- Major respiratory organizations recommend targeted genetic screening for alpha-1 in all patients with COPD⁵⁻⁹
- **AlphaID and AlphaID CONFIRM kits contain everything you need to screen for alpha-1**



Grifols is Committed to Supporting Patients with Alpha-1

Grifols is dedicated to advancing the understanding of alpha-1. Grifols has been the leader in alpha-1 screening for more than 30 years and continues to develop innovative ways to screen patients.

To begin screening with **AlphaID**, call your sales representative or visit **AlphaID.com**



References: **1.** Campos MA, Wanner A, Zhang G, Sandhaus RA. Trends in the diagnosis of symptomatic patients with α 1-antitrypsin deficiency between 1968 and 2003. *Chest*. 2005;128(3):1179-1186. **2.** Brantly ML, Paul LD, Miller BH, Falk RT, Wu M, Crystal RG. Clinical features and history of the destructive lung disease associated with alpha-1-antitrypsin deficiency of adults with pulmonary symptoms. *Am Rev Respir Dis*. 1988;138(2):327-336. **3.** Data on file, PROLASTIN DIRECT Program, Grifols. **4.** DeMeo DL, Sandhaus RA, Barker AF, et al. Determinants of airflow obstruction in severe alpha-1-antitrypsin deficiency. *Thorax*. 2007;62(9):806-813. **5.** American Thoracic Society/European Respiratory Society. American Thoracic Society/European Respiratory Society statement: Standards for the diagnosis and management of individuals with alpha-1 antitrypsin deficiency. *Am J Respir Crit Care Med*. 2003;168(7):818-900. **6.** CHEST Foundation. Learn About Alpha-1 Antitrypsin Deficiency. <https://foundation.chestnet.org/patient-education-resources/alpha-1-antitrypsin-deficiency/>. Last reviewed January 2019. Accessed April 10, 2023. **7.** Sandhaus RA, Turino G, Brantly ML, et al. The diagnosis and management of alpha1-antitrypsin deficiency in the adult. *Chronic Obstr Pulm Dis*. 2016;3(3):668-682. **8.** Global Initiative for Chronic Obstructive Lung Disease. *Pocket Guide to COPD Diagnosis, Management, and Prevention: 2020 Report*. 2022:1-53. **9.** World Health Organization. Global surveillance, prevention and control of chronic respiratory diseases: a comprehensive approach. 2017. https://www.who.int/gard/publications/GARD_Manual/en/. Accessed April 10, 2023. **10.** Carpenter MJ, Strange C, Jones Y, et al. Does genetic testing result in behavioral health change? Changes in smoking behavior following testing for alpha-1 antitrypsin deficiency. *Ann Behav Med*. 2007;33(1):22-28. **11.** Data on file, Alpha-1 Genetics Laboratory, Grifols. **12.** A1AT Genotyping Test Package Insert. Progenika Biopharma.